

Overview

The MTT690 is an automated tensile tester designed to measure the tensile properties of single hair fibres by stretching them to a specified percentage or to failure. The MTT690 was developed to overcome the low productivity associated with manual testing of mechanical properties of single fibres. The system is based on a circular sample cassette, which allows the automatic measurement up to 100 premounted fibre samples.

Principal benefits:

- Automated operations and analysis
- High throughput: 100 fibre rotary cassette
- Can incorporate additional module for wet testing
- Multi-tasking allows for simultaneous measurement, reducing resting time

Applications and claims:

- Strength claims
- Hydration claims
- Damage repair/alleviation claims




Metrology principle —

Hair fibre samples are mounted using brass crimps and placed onto a 100-slot rotary cassette, from which a pneumatically operated sample gripper picks up the sample. The gripper is mounted onto a load cell which measures the force being applied to the sample, which can be stretched to either a specific percentage or to failure. The MTT690 can also be integrated onto the automated sample loading system (ALS1500) with the FDAS770 to provide both dimensional and tensile data from a single fully automated measurement operation. Normalising tensile data with fibre cross-sectional area reduces data variability up to 80%, helping with group significant discrimination


Dedicated software - UvWin

The MTT690 is controlled using Dia-Stron UvWin software. UvWin offers several analysis options, e.g. 1 or 3 phase tensile analysis, hysteresis analysis and stress relaxation, as well as a number of integrated data processing tools. Raw data can also be exported as a text file, for use in Excel or other statistical packages. UvWin is compatible with the latest versions of the Windows OS.


Sample data and analysis—



Stress-strain curve created in UvWin software



Boxplot of the break stress, where the circle denotes the mean break stress



References —

Publications:


T.A. Evans, J.M. Altland, R. Burch and J. Horsager (2013): "Glycolic Acid No Longer Just for Skin – Changing the Internal Properties of Hair"; Cosmetics & Toiletries, Vol. 128, No. 11, November 2013, pages 816-820.

A. Dussaud, R. Bhavna, and L.T. Hui (2013): "Progressive hair straightening using an automated flat iron: Function of silicones"; Journal of Cosmetic Science, 64, 119–131.

Dan Istrate, Crisan Popescu, Meriem Er Rafik, Martin Möller (2013): "The effect of pH on the thermal stability of fibrous hard alpha-keratins"; Polymer Degradation and Stability 98 (2013) pages 542-549.

P. Morganti, M. Palombo, A. Cardillo, P. Del Ciotto, G. Morganti and G. Gazzaniga (2012): "Antidandruff and anti-oily Efficacy of Hair formulations with a Repairing and Restructuring activity. The Positive Influence of the Zn-Chitin Nanofibrils Complexes"; Journal of Applied Cosmetology, 30, 149-159 (July/December 2012).

M.G. Davis, J.H. Thomas, S. van de Velde, Y. Boissy, T.L. Dawson Jr, R. Iveson and K. Sutton (2011): "A novel cosmetic approach to treat thinning hair"; British Journal of Dermatology, 2011 Dec; 165 Suppl 3:24-30.



References —

Books:

Trefor Evans and R. Randall Wickett (2012): "Practical Modern Hair Science"; Chapter "Mechanical Properties of Hair", page 229; ISBN: 978-1932633931

Clarence R. Robbins (2012): "Chemical and Physical Behavior of Human Hair"; Chapter "The Physical Properties of Hair Fibers", page 537; 5th ed. 2012, XXIII, ISBN: 978-3-642-25610-3

Randy Schueller, Perry Romanowski (2011): "Hair Care Formulation for All Hair Types"; Chapter "Testing Polysilicone-19 for Hair Conditioning and UV Protection Claims", page 113; ISBN: 978-1932633863

Ulrike Blume-Peytavi, David A. Whitting (2008): "Hair Growth and Disorders"; Chapter "Mechanical Test of Hair Quality (Elasticity, Strength, Fragility)", page 142; ISBN: 978-3-540-46908-7

Examples of use in patent claims:

WO2018178056A1 Hair shaping material (Croda/University of Leeds) March 2018 US20140190507 Compositions And Methods For Enhancing The Structure Of Hair Fibers (Estee Lauder) July 2014

US8697040 Hair treatment compositions (Access Business Group) Apr 2014

WO2014025686 Compositions and Methods for Treating a Keratin Based Substrate (Ashland) Feb 2014 US8273332 Hair care product containing acetylpyridinium salts (Henkel) Sept 2012

US20120201776 Hair Relaxer (BASF) Aug 2012

US8221731 Continuous moisturization compositions (Aveda) July 2012

 $US8158116\,Method\,for\,treating\,hair\,damaged\,by\,color\,treatments\,(Ecolab/Lubrizol)\,April\,2012$

MTT690 Tensile Tester		
Extension range	0-120mm	
Speed range	1-2000mm/min	
Force range	0 to 20N (2000gmf)	
Force resolution	0.05gmf	
Displacement resolution	10µm	
Sample size	30mm	

Programmable Feature	es
----------------------	----

Methods

Strain/stress method with break detection Stress relaxation Hysteresis Creep

Content

PU1100 MTT690 Mechanical Unit UvWin software for Windows OS

UV1000

Power supply

Universal 85-265V AC 47-63Hz, 50W

Compressed air

Dry, clean compressed air. 4.5bar min, 20I/min

Computer Windows OS: 7 and 10, 1 x USB port

