

# Recommend the Suitable Sunscreen Factor



## Sun Protection Diagnostic SP 37



### pigmentation

- ideal device to analyse natural protection time (time somebody can expose to the sun unprotected) and sun screen needed according to sun intensity (location, season, reflection) and exposure time by pigmentation measurements of skin
- results are shown on a display and diode chains

The SP 37 is light, portable, battery operated, flexible and easy to handle.

It can be used easily on the counter, in a beauty or hairdressers salon, at the dermatologist or in the pharmacy to attract customers and to increase the sales of cosmetic products.

Please ask for our special transportation accessories.



	<1,000 m	>1,000 m
1	I	IV V
2	II III	V VI
3	IV V	VII VII
4	VIII VIII	VIII VIII



	<1,000 m	>1,000 m
1	II III	V VII
2	III IV	VI VII
3	V VI	VII VII
4	VIII VIII	VIII VIII



	<1,000 m	>1,000 m
1	IV V	VII VII
2	V VI	VII VIII
3	VII VII	VIII VIII
4	VIII VIII	VIII VIII

Technical Data:  
 Power Supply: battery operated (9 V), optional external power supply (100-240VAC/9 VDC) Dimensions: 17 x 20 x 5.5 cm; Weight: 0.5 kg;  
 Measurement Principle: melanin by absorption/reflection

Technical changes may be made without prior notice.